

ZEN (blue edition) 3.1

Image Analysis



Dr. Marion Lang
Product Management
2019-11-04

Getting Started

Simulate Hardware




Copy **Active Configuration.xml** and **CZIS_Cameras.xml** to
C:\ProgramData\Carl Zeiss\MTB2011\2.16.0.9


Start MTB → Check “Simulate”

C:\Program Files\Carl Zeiss\MTB 2011 - 2.11.0.0\MTB Configuration\MTBConfig.exe


Check that there is a demo camera




Add Python to the PATH variable




Alternatively: In the script specify the whole path to the python.exe



Registry Entry



Simulate the Acquisition




Set path to demo files, e.g.
...\CCD_Testimages_Folder\96_well

Familiarize with ZEN Application and Image Analysis



Task: Acquire some images and set up a simple image analysis with cell counting.


Testdata: ...\\Simple Cell Count\\Cell Count.czi




Familiarize with ZEN Application and Image Processing



Task: Apply an Image Processing Function to the data





OAD Macro Development: Getting Help

Macro object model



Macro Object Model

The screenshot shows the Macro Editor interface with the 'Help' menu open, highlighting 'Macro Object Model...'. The Macro Editor window displays the code `img = Zen.Application.ActiveDocument`. The Macro Object Model Documentation window is open, showing the 'Segmentation.Threshold Method' page. The page includes a search bar with 'threshold' entered, a list of search results, and an 'Overload List' table.


Name	Description
<code>Threshold(ZenImage)</code>	This function performs a segmentation based on the definition of a brightness range for the regions to be segmented.
<code>Threshold(ZenImage, Boolean)</code>	This function performs a segmentation based on the definition of a brightness range for the regions to be segmented.
<code>Threshold(ZenImage, ThresholdSetting)</code>	Call of function with ThresholdSetting.

OAD Macro Development: Getting Help

Intellisense Auto-Completion



Intellisense Auto-Completion
type **Zen / Zen.**



OAD Macro Development: Getting Help

OAD Forum



<http://forums.zeiss.com/microscopy/community/>

Open Application Development (OAD) for ZEN (blue edition)

Threads / Posts Last Post

You have a special application that demands functionality beyond ZEN? Use the integrated OAD (Open Application Development) environment of ZEN (blue edition). With OAD you create your own macro solution based on the well-established Python language. Benefit from the simple access to a vital set of ZEN functions and the ability to include libraries such as the .Net Framework. Join the OAD community to discuss macros and help other users to develop their ultimate solution.



Instrument Control (2 Viewing)

Discuss macros to control the hardware functions of your ZEISS microscope

Threads: 55
Posts: 201

SmartProof5 - Control... »
by Carl Zeiss Microscopy 3
10-24-2019, 08:59 AM



Image Acquisition (1 Viewing)

Post your acquisition-related questions and macros here

Threads: 42
Posts: 113

SmartProof5 - Import Tiles... »
by Carl Zeiss Microscopy 3
10-24-2019, 09:19 AM



Image Handling (1 Viewing)

Your place to discuss general handling of images and the CZI file format

Threads: 53
Posts: 146

Problems using czi image... »
by Fredrik
10-31-2019, 10:09 AM



Image Processing

Discuss questions and projects related to processing of imaging data here

Threads: 56
Posts: 189

ZeissImageLib and .Net Core »
by MosGeo
10-31-2019, 01:59 PM



Measurement and Analysis (1 Viewing)

Share your programming ideas for various measurement and analysis tasks here

Threads: 24
Posts: 66

Extract rotated rectangle... »
by Fredrik
10-22-2019, 03:54 PM

OAD Macro Development: Getting Help

Sample Macros on Github



<https://github.com/zeiss-microscopy/OAD>

Collection of tools and scripts useful to automate microscopy workflows in ZEN Blue using Python and Open Application Development tools.

python microscopy scripting automation imaging machine-learning open-source image-analysis zen zen-blue acquisition
python-script workflow oad image-processing

145 commits

1 branch

0 releases

2 contributors

GPL-3.0

Branch: master ▾

New pull request

Find file

Clone or download ▾

sebi06 updates










Latest commit 405a983 4 days ago

OAD Macro Development: Getting Help




Sample Macros on ZEN DVD



Manuals → ZEN (blue edition) → OAD content → Experiment Feedback

 Exp_Scripts_only Type: Folder	Date modified: 30.10.2019 16:31
 Exp_with_Scripts Type: Folder	Date modified: 30.10.2019 16:31
 Feedback_Script_Reader Type: Folder	Date modified: 30.10.2019 16:31
 Fiji Type: Folder	Date modified: 30.10.2019 16:31
 Image_Analysis Type: Folder	Date modified: 30.10.2019 16:31
 Python_Scripts Type: Folder	Date modified: 30.10.2019 16:31
 SoundFiles Type: Folder	Date modified: 30.10.2019 16:31
 Testimages Type: Folder	Date modified: 30.10.2019 16:31
 Experiment_Feedback_Tutorial_3.1.pdf Type: Adobe Acrobat Document	Date modified: 28.10.2019 08:55 Size: 5,53 MB → 5,25 MB


Other → OAD Examples

 OAD Macro Collection Type: Folder	Date modified: 30.10.2019 16:32
 OAD Scripts and Tools Type: Folder	Date modified: 30.10.2019 16:32
 Disclaimer.txt Type: Text Document	Date modified: 28.10.2019 08:54 Size: 634 bytes → 297 bytes

Finding ZEN Folders



Help → About ZEN → Show ZEN Information → Folders → Doubleclick e.g. on UserDocuments



2_5d render series	05.10.2017 09:58	Dateiordner
2_5d render settings	05.10.2017 09:58	Dateiordner
3Dxl render series	05.10.2017 09:58	Dateiordner
3Dxl render settings	29.03.2019 10:03	Dateiordner
Adaptive Focus Point Settings	29.03.2019 10:03	Dateiordner
Automated Photomanipulation Settings	22.08.2019 13:51	Dateiordner
Available Features	29.03.2019 10:03	Dateiordner
Available Graphic Elements	29.03.2019 10:03	Dateiordner
Camera Settings	29.03.2019 10:03	Dateiordner
Chart Settings	05.10.2017 09:58	Dateiordner
Color Schemes	05.10.2017 09:58	Dateiordner
Display Settings	29.03.2019 10:03	Dateiordner
Experiment Setups	02.11.2019 13:58	Dateiordner
Graphics	29.03.2019 10:03	Dateiordner
Guided Acquisition Settings	01.10.2019 16:16	Dateiordner
Hardware Settings	29.03.2019 10:03	Dateiordner
Image Analysis Settings	31.10.2019 19:29	Dateiordner
ImageJ	29.03.2019 10:03	Dateiordner
Interactive Measurement Sequences	29.03.2019 10:03	Dateiordner
Live Scan Profiles	29.03.2019 10:03	Dateiordner
Live Scan Sample Holder Templates	25.10.2017 11:59	Dateiordner
Macros	02.11.2019 14:58	Dateiordner
Measure Features Selection	29.03.2019 10:03	Dateiordner
Measure Features Subset	29.03.2019 10:03	Dateiordner
Model-Repository	25.07.2019 10:48	Dateiordner
Multicaption Blocking Settings	21.09.2019 09:55	Dateiordner


OAD Macro Development: Getting Help

Macro recorder







Note: Not everything can be recorded!

But that does not mean, there is not command for a specific functionality






Almost everything is done via settings. E.g. Experiment



Setting




XML-File




Acquisition Experiment Manager

Almost everything is done via settings. E.g. Processing



Setting




XML-File




Processing

Almost everything is done via settings. E.g. Measurement



Setting





Processing

XML-File

Almost everything is done via settings. E.g. Image Analysis



Analysis



Setting



XML-File

Hardware:


Zen.Devices.

Lists all devices/methods that are directly accessible, e.g.


posX = Zen.Devices.Stage.ActualPositionX

mag = Zen.Devices.ObjectiveChanger.Magnification

Zen.Devices.



Hardware Setting:



Read & Save Hardware Setting



#Read hardware setting:

```
HWS = Zen.Devices.ReadHardwareSetting()
```

#Save hardware setting:

```
HWS.SaveAs(filename)
```

#Create an empty container:

```
HWS = ZenHardwareSetting()
```

#Load an existing Hardwaresetting:

```
HWS.Load('filename')
```

Default file path:

```
C:\Users\...\Documents\Carl
```

```
Zeiss\ZEN\Documents\Hardware Settings
```

→ OAD_Training_Get HWS

Get Components from the Hardware Setting



```
# Get the IDs of every component in the Hardware Setting:
HWSIDs = HWS.GetAllComponentIds()

# go through all the Component IDs
for ID in HWSIDs:
    print(ID)

# for each Component ID get the Parameter names
HWSParams = HWS.GetAllParameterNames(ID)

# go through all Parameters
for HWSParam in HWSParams:
    print("\t" + HWSParam)
```

```
AcquisitionFrame.Max
LiveImageFrame.Max
Frame.Enabled
CameraIdentifier
CameraIdentifier.Type
CameraIdentifier.Default
CameraIdentifier.GuiHint
CameraIdentifier.DisplayOrder
CameraIdentifier.Label
TheoreticalTotalMagnification
TotalMagnification
DefaultScalingUnit
RoiCenterOffsetX
RoiCenterOffsetY
TotalAperture
SoftwareAutofocus
    IsInitialized
    IsRunning
    AutofocusResult
    AutofocusSignalQuality
    AutofocusAnalysisData
Autofocus
    IsInitialized
    IsRunning
    AutofocusResult
    AutofocusSignalQuality
    AutofocusStrategyMode
    AutofocusSignalQualityThreshold
ImmersionAdapter
    IsAutoImmersionEnabled
    IsImmersionEnabled
    IsImmersionAvailable
    IsImmersionPrimed
```


Set a parameter and apply



#Get a parameter:

```
stage_speed = HWS.GetParameter('componentId',  
    'parameterName')
```

#Set a parameter:

```
HWS.SetParameter('componentId', 'parameterName',  
    parameterValue)
```

#Apply the setting:

```
Zen.Devices.ApplyHardwareSetting(HWS)
```

#Apply a single value:

```
HWS = ZenHardwareSetting()  
HWS.SetParameter('set parameters')  
HWS.SetActive()  
HWS.SaveAs('filename')
```

- Start macro recorder.
- Load an Image
- Use IP „Threshold“

```
# ***** Recorded Code Block *****  
image1 = Zen.Application.LoadImage("C:\\testdata\\Koala_Graphics.czi", True)  
Zen.Application.Documents.Add(image1)  
image2 = Zen.Processing.Segmentation.Threshold(image1, 53, 181, 198, 127, 0,  
199, True, False, False, True)  
Zen.Application.Documents.Add(image2)  
# ***** End of Code Block *****
```



- Create a setting for the IP „Threshold“
- Apply the IP with threshold to the image
- Add the image to the document collection

Get the image

```
image1 = Zen.Application.Documents.GetByName("Koala_Graphics.czi")
```

Create a setting object

```
set = Zen.Processing.Segmentation.Settings.ThresholdSetting()
```

Load the setting


```
set.Load("C:\\Users\\M1MALANG\\Documents\\Carl  
Zeiss\\ZEN\\Documents\\Processing  
Settings\\Threshold\\koala_threshold.czip")
```

Run IP

```
image2 = Zen.Processing.Segmentation.Threshold(image1, set, True)
```

Add result to image documents


```
Zen.Application.Documents.Add(image2)
```



Allows to run a macro before or after experiment execution
Will be part of the experiment xml.

- Create a macro, e.g. „invert“
- Run it automatically at the end of the experiment

```
# invert current image  
image1 = Zen.Application.ActiveDocument  
image2 = Zen.Processing.Binary.Invert(image1)  
Zen.Application.Documents.Add(image2)  
image1.Close()
```



OAD

Customize Application



Tools → Customize Application




Image Analysis with OAD

Run an image analysis, generate tables



```
img = Zen.Application.ActiveDocument
ias = ZenImageAnalysisSetting()
ias.Load("cell_count")
```

```
# run image analysis
```

```
img = Zen.Analyzing.Analyze(img, ias)
```


„Analyze“ stores the results in the czi file

```
# create result table
```

```
table = Zen.Analyzing.CreateRegionTable(img)
table.Save("C:\\temp\\cell_count_region.csv")
Zen.Application.Documents.Add(table)
```

```
img = Zen.Application.ActiveDocument
ias = ZenImageAnalysisSetting()
ias.Load("cell_count")

#create labeled images
img_out = ZenImage()
Zen.Analyzing.AnalyzeToImage(
img, img_out, ias,
ZenAnalyzerLabelImageMode.RegionUniqueColor,
ZenAnalyzerLabelImageInitMode.CopyAllChannels,
Zen.PixelType.Bgr24)
Zen.Application.Documents.Add(img_out)
```



Start External Programs

Start Fiji and load result in ZEN



```
img = ZenImage()  
img.Load("C:\\Training\\CZI_DimorderTZC.czi")  
  
filename = img.FileName  
exeloc = 'C:\\Fiji_Sebi\\ImageJ-win64.exe'  
  
# specify the Fiji macro that will be applied  
macro = r'-macro C:\\Training\\Open_CZI_and_MaxInt_Save.ijm'  
  
# define parameters for the Fiji macro  
params = macro + ' ' + filename  
  
# start Fiji, open the data set and execute the macro  
app = Process();  
app.StartInfo.FileName = exeloc  
app.StartInfo.Arguments = params  
app.Start()  
app.WaitForExit()
```

→ OAD_Training_Fiji

Start External Programs

Start Fiji and load result in ZEN



ImageJ Macro

Open_CZI_and_MaxInt_Save.ijm

```
1
2 name = getArgument;
3 if (name=="") exit("No argument!");
4
5 run("Bio-Formats Importer", "open=[" + name + "] autoscale color_mode=Default open_all_series view=Hyperstack stack_order=XYCZT");
6 // get dimension from original image
7 getDimensions(width,height,channels,slices,frames);
8 // apply additional operation
9 run("Z Project...", "start=1 stop=" + slices + " projection=[Max Intensity]"); // do maximum intensity projection
10 run("Fire"); // apply special LUT
11 saveAs("PNG", "C:/Training/fiji.png");
12 run("Quit");
```

```
savename = "C:\\Training\\fiji.png"
```

```
if File.Exists(savename):
    fiji_result = Zen.Application.LoadImage(savename, False)
    Zen.Application.Documents.Add(fiji_result)
else:
    print 'Saved figure not found.'
```

Application Example: Translocation


Inhibition of nuclear export of Forkhead



Unhibited export

Add Wortmannin / LY294002

Inhibited export



OAD Script for automatization

Start Python Script for Calculation and Plotting



Tasks to perform:



1. Load load the image file (*.csv) and image anlaysis setting (*.czias)
2. Run the image analysis
3. Extract the image analysis results as *.csv
4. Start the python script (test_wellplate_from_ZEN.PY)



5. Read in data
6. Calculate the translocation Ratio
7. Generate heatmaps for different features (e.g. Translocation Ratio)
8. Save heatmaps as PNG files



9. Load PNG files in ZEN

Testdata: Translocation_comb_96_5ms.czi

Image Analysis Setting: Translocation_26.czias

```
from System.Diagnostics import Process
from System.IO import File, Path, Directory
import time
```

→ Translocation Plot Pandas

```
# define the external plot script or tool
pythonexe = r'C:\Program Files\Carl Zeiss\ZeissPython\Py20190211\env\python.exe' # Zeisspython
# requires progressbar2 library: "C:\Program Files\Carl Zeiss\ZeissPython\Py20190211\env\python.exe" -m pip install
progressbar2
script = r'C:\...\test_wellplate_from_ZEN.PY'


# load image and add it to ZEN
image_to_analyze = r'C:\...\Translocation_comb_96_5ms.czi'
image = Zen.Application.LoadImage(image_to_analyze)
Zen.Application.Documents.Add(image)

# get the image path
outputpath = Path.GetDirectoryName(image_to_analyze)
resultname = Path.GetFileNameWithoutExtension(image.Name)
```

Features defined in the image analysis setting:

ID parent | ID | Image Scene container Name | Image Scene Row | Image Scene Column | Index # | Area | NucMeanDapi
| NucMeanGFP | RingMeanGFP | RingArea

```
# Load image analysis setting and perform image analysis
iasfilename = r'C:\...\Image Analysis Settings\Translocation_2.czias'
ias = ZenImageAnalysisSetting()
ias.Load(iasfilename)
Zen.Analyzing.Analyze(image,ias)
```



```
# For ZOI-Image Analysis Settings need to get the results for the Primary Objects
# Create data list with results for each primary object
table_single = Zen.Analyzing.CreateRegionTable(image, "Primary Object")
#Zen.Application.Documents.Add(table_single)

# Save both data list as CSV files
table_single_filename = Path.Combine(outputpath, resultname + '_Single.csv')
table_single.Save(table_single_filename)


# close the image and image analysis setting
image.Close()
ias.Close()
```

	A	B	C	D	E	F	G	H	I	J	K	L
1	ParentID::ID	ID::ID!!!	ImageScene	ImageScene	ImageScene	Index::Index	Area::Area!!	IntensityMe	IntensityMe	CopyRingInt	CopyRingArea::Ring Area!!R	
2							µm²	Gray	Gray	Gray	µm²	
3	269	270 A1		1	1	1	2.64	133.337121	88.0530303	141.167832	1.43	
4	271	272 A1		1	1	2	1.84	79.2880435	0	0.02962963	1.35	
5	273	274 A1		1	1	3	3.04	168.177632	9.81907895	9.22058824	1.36	
6	275	276 A1		1	1	4	2.81	106.231317	38.7935943	46.6071429	1.96	
7	277	278 A1		1	1	5	2.47	146.546559	52.7773279	65.0347222	1.44	
8	279	280 A1		1	1	6	2.01	96.6865672	33.4079602	68.3777778	1.35	
9	281	282 A1		1	1	7	1.58	94.3734177	18.7468354	28.1842105	1.14	
10	283	284 A1		1	1	8	1.81	87.6574586	31.0607735	53.5683453	1.39	
11	285	286 A1		1	1	9	2.25	114.826667	0	0	1.36	
12	287	288 A1		1	1	10	2.72	95.4669118	0.11397059	0.625	1.44	
13	289	290 A1		1	1	11	1.95	96.4564103	13.2974359	22.0735294	1.36	
14	291	292 A1		1	1	12	2.23	104.753363	0.07623318	0.10218978	1.37	

```
# define the actual CSV file and the parameters
csvfile = Path.Combine(outputpath, table_single_filename)

# this depends on the actual CZIAS and the import of the CSV table in python
#parameter2display = 'CellCount'
parameter2display = 'Ratio'
params = ' -f ' + csvfile + ' -w 96' + ' -p ' + parameter2display + ' -sp False -dpi 100 -xlsx
True'

# start the data display script as an external application
Process.Start(script, params)
```




```
C:\Anaconda3\python.exe
CSV Filename: C:\testdata\Broadinstitute\Translocation_conf_96_Sms_Single.csv
PlateType: 96
Parameter to display: Ratio
DPI: 100
Columns : Index<['ParentID::ID of the parent!!I', 'ID::ID!!I',
'ImageSceneContainerName::Image Scene Container Name ',
'ImageSceneRow::Image Scene Row Index!!I',
'ImageSceneColumn::Image Scene Column Index!!I', 'Index::Index!!I',
'Area::Area!!R',
'IntensityMean_DAPI::Intensity Mean Value of channel 'DAPI'!!R',
'IntensityMean_EGFP::Intensity Mean Value of channel 'EGFP'!!R',
'CopyRingIntensityMean1::Ring Mean Intensity 1!!R',
'CopyRingArea::Ring Area!!R' ],
dtype='object')
Number of Object Parameters: 6
-----
wellID_key : WellID
Found keys:
Index<['ParentID', 'ID', 'WellID', 'RowID', 'ColumnID', 'Index', 'Area',
'NucMeanDapi', 'NucMeanGFP', 'RingMeanGFP', 'RingArea', 'Ratio' ],
dtype='object')
-----
FloatProgress(value=1.0, description='Processing Wells', max=96.0, min=1.0)
```

```
# define filenames of PNG files
savename_all = Path.Combine(Path.GetDirectoryName(image_to_analyze),
Path.GetFileNameWithoutExtension(image_to_analyze) + '_Single_HM_all.png')
savename_single = Path.Combine(Path.GetDirectoryName(image_to_analyze),
Path.GetFileNameWithoutExtension(image_to_analyze) + '_Single_HM_' + parameter2display + '.png')
print 'Showing saved figure in ZEN.'

# Check if filename exists, and Load images in ZEN
if File.Exists(savename_all):
    plotfigure1 = Zen.Application.LoadImage(savename_all, False)
    plotfigure2 = Zen.Application.LoadImage(savename_single, False)
    Zen.Application.Documents.Add(plotfigure1)
    Zen.Application.Documents.Add(plotfigure2)
else:
    print 'Saved figure not found.'

print 'Done.'
```



macrolib_2.py

```
def greet_user(username):  
    print("Hello, " + username + "!")
```

Macro:

```
import macrolib_2  
Zen.Application.MacroEditor.ClearMessages()  
  
macrolib_2.greet_user('Marion')
```

